

April 27th, 2015

Now Introducing, Medgyn's Endosampler II™!

Same technology as the Endosampler , only better! Our new Endosampler II™ is designed with two suction holes at tip for larger sampling size of tissue and reduce the force of suction for better comfort.

This device is designed to collect clearly-differentiated endometrial tissue without anesthesia, in an outpatient setting or as an office procedure. This device improves efficiency, convenience and reliability of the endometrial tissue collection and to reduce the physical discomfort of the patient. It provides liberal histologic specimens, surgical ease and is procedurally cost-effective for routine use. MedGyn's Endosampler II™ is recommended for both pre- and post-menopausal women for diagnostic purposes.

MedGyn's Endosampler II™ is available in a semi-rigid 3mm curette with a safe luer lock 10 ml syringe.